

D3-Triclosan (2,4-dichlorophenoxy D3) 100 µg/mL in Acetonitrile

Art. ID HPC-675735-1ML
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|----------------|---------------------|-------|--------|--------|
| | Triclosan D3 (2,4-dichlorophenoxy D3) | [1020719-98-5] | 100 | µg/mL | | |