

D3-Mecoprop 10 µg/mL in Acetone

Art. ID HPC-677823-10ML
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|---------------|---------------------|-------|--------|--------|
| | Mecoprop D3 (ring) | [352431-15-3] | 10 | µg/mL | | |