

Custom Solution - VOC - Mix 11 in Methanol

Art. ID LM24-CUS-10103
Unit 5 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 18 months (originally closed bottle)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,2,3,3-Tetramethylsucc inonitrile | [3333-52-6] | 20 | µg/mL | | |
| | Aniline | [62-53-3] | 20 | µg/mL | | |
| | 4-Fluoro-N-isopropylani line | [70441-63-3] | 20 | µg/mL | | |
| | 3-Cyanopyridine | [100-54-9] | 20 | µg/mL | | |
| | 2-Chloroaniline | [95-51-2] | 20 | µg/mL | | |
| | 4-Chloroaniline | [106-47-8] | 20 | µg/mL | | |
| | 2-Chloro-5-methylpyridi ne | [18368-64-4] | 20 | µg/mL | | |
| | 2,3-Dichloro-1-propanol | [616-23-9] | 20 | µg/mL | | |
| | 1,4-Thioxane | [15980-15-1] | 20 | µg/mL | | |