

## **Custom Solution - VOC - Mix 11 in Methanol**

Art. ID LM24-CUS-10105  
Unit 5 x 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

### Description

Conservation 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2,3,3-Tetramethylsucc inonitrile	[3333-52-6]	100	µg/mL		
	Aniline	[62-53-3]	100	µg/mL		
	4-Fluoro-N-isopropylani line	[70441-63-3]	100	µg/mL		
	3-Cyanopyridine	[100-54-9]	100	µg/mL		
	2-Chloroaniline	[95-51-2]	100	µg/mL		
	4-Chloroaniline	[106-47-8]	100	µg/mL		
	2-Chloro-5-methylpyridi ne	[18368-64-4]	100	µg/mL		
	2,3-Dichloro-1-propanol	[616-23-9]	100	µg/mL		
	1,4-Thioxane	[15980-15-1]	100	µg/mL		