

Nanosphere Size Standard (Nominal: 400 nm) (HS-Series), LOT: HS393.115 (NIST traceable)

Art. ID BSP-HS0040-20
Unit 20 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------|------------|---------------------|------|--------|--------|
| | Certified Mean Diameter | | 401 | nm | | |
| | Size CV | | 4,7 | % | | |
| | Solids | | 2 | % | | |