

PharamCount(TM) TriCal - Size standard and count control for 3-point calibration of autodiluters (Set of 6 particle suspensions and certificate)

Art. ID BSP-20825
Unit 25 mL each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Specifications: 0,8 µm ± 5 %: nom. 1.200.000 P article/ mL, 10 µm ± 5 %: nom. 800.000 Particl e/ mL, 5 µm ± 5 %: nom. 400.000 Particle/ mL						