

Zeta potential control - Suspension of modified polystyrene for reproducibly measuring the zeta potential - Negatively charged, 0,02% solids, in a ready-to-use syringe with Luer-lock cap (18 syringes)

Art. ID BSP-70318C
Unit 3,3 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------|------------|---------------------|------|--------|--------|
| | Zeta potential (nominal) | | (- 50) | mV | | |
| | Nominla mean size | | 200 | nm | | |