

Zeta potential control - Suspension of modified polystyrene for reproducibly measuring the zeta potential - Negatively charged, 0,02% solids, in a ready-to-use syringe with Luer-lock cap (18 syringes)

Art. ID BSP-70318C
Unit 3,3 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zeta potential (nominal)		(- 50)	mV		
	Nominla mean size		200	nm		