

Custom Solution - Lot A 5000 µg/mL each of Mineral oil type A and Mineral oil type B in n-Heptane

Art. ID LM24-CUS-70443
Unit 5x10 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1206 II E2

Description

Conservation 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Lot A	Mineral oil type A	[n/a]	5000	µg/mL		
Lot A	Mineral oil type B	[n/a]	5000	µg/mL		