

Custom Solution - VOC Mix 4 in Methanol

Art. ID LM24-CUS-10113

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Ethoxyethanol	[110-80-5]	100	µg/mL		
	Tetrachloroethene	[127-18-4]	100	µg/mL		
	2-Methylphenol (o-Creso l)	[95-48-7]	100	µg/mL		
	Ethylbenzene	[100-41-4]	100	µg/mL		