

### **Custom Solution - VOC - Mix 4 500 µg/ml in Methanol**

Art. ID LM24-CUS-10109  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### **Description**

Conservation 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1-Methyl-2-pyrrolidon	[872-50-4]	500	µg/mL		
	N,N-Dimethylacetamide	[127-19-5]	500	µg/mL		
	N,N-Dimethylformamide	[68-12-2]	500	µg/mL		
	Formamide	[75-12-7]	500	µg/mL		