

### **Custom Solution - VOC - Mix 3 je 100 µg/ml in Methanol**

Art. ID LM24-CUS-10111  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tris(2-chloroethyl)phosphate	[115-96-8]	100	µg/mL		
	Tetrabromobisphenol A	[79-94-7]	100	µg/mL		
	2,2-Bis-(bromomethyl)-1,3-propanediol	[3296-90-0]	100	µg/mL		