

Leaded Brass, Chill Cast, Certified Reference Material (CRM), disc 40 mm x 15 mm

Art. ID MBH-31X7835.11A
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	0,1522	%		
	Lead (Pb)	[7439-92-1]	1,695	%		
	Iron (Fe)	[7439-89-6]	0,192	%		
	Nickel (Ni)	[7440-02-0]	0,1007	%		
	Aluminium (Al)	[7429-90-5]	0,908	%		
	Copper (Cu)	[7440-50-8]	65,75	%		
	Manganese (Mn)	[7439-96-5]	0,0101	%		
	Antimony (Sb)	[7440-36-0]	0,0122	%		
	Cobalt (Co)	[7440-48-4]	0,0065	%		
	Arsenic (As)	[7440-38-2]	0,143	%		
	Phosphorus (P)	[7723-14-0]	0,024	%		
	Silicon (Si)	[7440-21-3]	0,01	%		
	Bismuth (Bi)	[7440-69-9]	0,011	%		
	Tellurium (Te)	[13494-80-9]	~0,009	%		
	Zinc (Zn)	[7440-66-6]	30,98	%		