

Waste water - Trace metals, LOW-AI 200 ng/mL

Art. ID SPS-WW1(LOW-AI 200 ng/mL)
Unit 6 x 50 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	200	ng/mL		
	Arsenic (As)	[7440-38-2]	100	ng/mL		
	Cadmium (Cd)	[7440-43-9]	20	ng/mL		
	Cobalt (Co)	[7440-48-4]	60	ng/mL		
	Chromium (Cr)	[7440-47-3]	200	ng/mL		
	Copper (Cu)	[7440-50-8]	400	ng/mL		
	Iron (Fe)	[7439-89-6]	1000	ng/mL		
	Manganese (Mn)	[7439-96-5]	400	ng/mL		
	Nickel (Ni)	[7440-02-0]	1000	ng/mL		
	Phosphorus (P)	[7723-14-0]	1000	ng/mL		
	Lead (Pb)	[7439-92-1]	100	ng/mL		
	Vanadium (V)	[7440-62-2]	100	ng/mL		
	Zinc (Zn)	[7440-66-6]	600	ng/mL		