

PAH Standard Solution - 16 components 2000 µg/mL each of in Methylene chloride (Dichloromethane)/Benzene (1/1)

Art. ID LM24-OP-4010095
Unit 1 mL
Deliverydetails De Minimis UN1992 (6.1) III E1

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene	[83-32-9]	2000	µg/mL		
	Acenaphthylene	[208-96-8]	2000	µg/mL		
	Anthracene	[120-12-7]	2000	µg/mL		
	Benzo(a)anthracene	[56-55-3]	2000	µg/mL		
	Benzo(a)pyrene	[50-32-8]	2000	µg/mL		
	Benzo(b)fluoranthene	[205-99-2]	2000	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	2000	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	2000	µg/mL		
	Chrysene	[218-01-9]	2000	µg/mL		
	Dibenz(a,h)anthracene	[53-70-3]	2000	µg/mL		
	Fluoranthene	[206-44-0]	2000	µg/mL		
	Fluorene	[86-73-7]	2000	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	2000	µg/mL		
	Naphthalene	[91-20-3]	2000	µg/mL		
	Phenanthrene	[85-01-8]	2000	µg/mL		
	Pyrene	[129-00-0]	2000	µg/mL		