

Custom Solution - 5 Components 50 mg/L each of in 2% HNO₃

Art. ID LM24-CUS-70463
Unit 2 x 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation: 12 months (originally closed bottle), 6 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Iron (Fe)	[7439-89-6]	50	mg/L		
	Manganese (Mn)	[7439-96-5]	50	mg/L		
	Barium (Ba)	[7440-39-3]	50	mg/L		
	Uranium (U)	[7440-61-1]	50	mg/L		