

Custom Solution - 1 Components Potassium (K) 50 g/L each of in 5% HNO₃

Art. ID LM24-CUS-70471
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]		50	g/L	