

### **Ni/Co, CRM, 1/4 section of 88 mm x 20 mm**

Art. ID IMZ-206  
Unit each  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 4,99                | %    |        |        |
|                  | Boron (B)         | [7440-42-8] | 0,001               | %    |        |        |
|                  | Carbon (C)        | [7440-44-0] | 0,035               | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 5,37                | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 9,78                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,0003              | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,036               | %    |        |        |
|                  | Hafnium (Hf)      | [7440-58-6] | 0,36                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,002               | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 0,064               | %    |        |        |
|                  | Niobium (Nb)      | [7440-03-1] | 0,023               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 62,24               | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,0007              | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,0004              | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,028               | %    |        |        |
|                  | Tantalum (Ta)     | [7440-25-7] | 11,95               | %    |        |        |
|                  | Titanium (Ti)     | [7440-32-6] | 1,36                | %    |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 4,08                | %    |        |        |
|                  | Zirconium (Zr)    | [7440-67-7] | 0,004               | %    |        |        |