

Custom Solution - 4 components in n-Hexane

Art. ID LM24-CUS-70476
Unit 1 x 1 mL
Deliverydetails De Minimis UN1208 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mineral oil type A	[n/a]	5	mg/L		
	Mineral oil type B	[n/a]	5	mg/L		
	n-Tetracontane	[4181-95-7]	2	mg/L		
	n-Decane	[124-18-5]	2	µL/L		