

Listeria monocytogenes dna agarose plug

Art. ID	ERM-AD624
Unit	vial
Deliverydetails	No Dangerous Good /not restricted

Description

Each unit consists of one vial containing one agarose plug for PFGE with embedded undigested DNA of Listeria monocytogenes. The plugs were produced and certified under the responsibility of the European Commission's Joint Research Centre. Listeria monocytogenes H2446 strain was used as the starting material for the production of the plugs. The plugs are stored in microvials containing TRIS-EDTA buffer solution

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please see certificate						