

Listeria monocytogenes dna agarose plug

Art. ID	ERM-AD624
Unit	vial
Deliverydetails	No Dangerous Good /not restricted

Description

Each unit consists of one vial containing one agarose plug for PFGE with embedded undigested DNA of *Listeria monocytogenes*. The plugs were produced and certified under the responsibility of the European Commission's Joint Research Centre. *Listeria monocytogenes* H2446 strain was used as the starting material for the production of the plugs. The plugs are stored in microvials containing TRIS-EDTA buffer solution

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please see certificate						