

## **Tool Steel - disc (38 mm x 25 mm)**

Art. ID ECRM-S 268-1-D  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

### **Description**

The ECRM-S 268-1-D is available as 38 mm diameter disc 25 mm thick.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Carbon (C)        | [7440-44-0] | 1,134 ± 0,007       | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,373 ± 0,007       | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,293 ± 0,005       | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,0209 ± 0,0006     | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,0154 ± 0,0004     | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 4,578 ± 0,021       | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 3,208 ± 0,019       | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,1437 ± 0,0022     | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,0062 ± 0,0003     | %    |        |        |
|                  | Boron (B)         | [7440-42-8] | 0,0009 ± 0,0002     | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 0,029 ± 0,0005      | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,1232 ± 0,002      | %    |        |        |
|                  | Nitrogen (N)      | [7727-37-9] | 2,03 ± 0,036        | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,0078 ± 0,0002     | %    |        |        |
|                  | Vanadium (V)      | [7440-62-2] | 8,478 ± 0,052       | %    |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 3,707 ± 0,029       | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,0017 ± 0,0002     | %    |        |        |