

## **Secondary Ferrite Number Standard - Low Range**

Art. ID	NIST-8480
Unit	10x12x20
Deliverydetails	No Dangerous Good /not restricted

### Description

This Reference Material (RM) is intended for the calibration of instruments used to measure weld metal ferrite content in accordance with American National Standards Institute (ANSI), American Welding Society (AWS) and International Standards Organization (ISO) normative standards ANSI/AWS A4.2 and ISO 8249. The ferrite content is measured in terms of a method defined quantity, the Ferrite Number (FN). NIST measurements were made in accordance with ANSI/AWS A4.2.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please aks for the current certificate						