

Microbial Genomic DNA Standards for Sequencing Performance Assessment

Art. ID	NIST-8375
Unit	4 vials, 1 each level
Deliverydetails	No Dangerous Good /not restricted

Description

This Reference Material (RM) is intended for validation, optimization, process evaluation, and performance assessment of microbial whole genome sequencing. Particular attention was devoted to assessing the genome sequence. A unit of RM 8375 consists of four vials. Each vial contains a different microbial genomic DNA sample (MG-001 Salmonella Typhimurium LT2, MG-002 Staphylococcus aureus, MG-003 Pseudomonas aeruginosa, and MG-004 Clostridium sporogenes).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please aks for the current certificate						