

### **Custom Solution - Stearyl stearate 1000 mg/L in Methanol**

Art. ID LM24-CUS-70488  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Stearyl stearate	[2778-96-3]	1000	mg/L		