

## **xrPress ET25 - MOTORIZED Hydraulic Presses**

Art. ID                      LM24-PRESS-ET25  
Unit                         68kg

### Description

The MotorT25, is a power assisted hydraulic press operating to 25 Tons. The press has been designed to handle a wide variety of pressing applications, including XRF and IR sample preparation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
<b>Key features</b>						
Large working distance between pressing surfaces						
Low noise operation						
Integral high clarity PETG safety Guards						
Max. Load on Ram (Piston) T25: 25 Tons						
Load Range (0.5Ton increments) T25: 2 – 25 Ton						
Top Bolster Diameter: 400mm						
Top Lead Screw Vertical travel: 90mm						
Ram (Piston) Bolster Diameter: 100mm						
Ram (Piston) Stroke: 25mm						
Max/Min Dist. between Pressing Faces: 165mm – 70mm						
<b>Technical specifications</b>						
Dimensions (H x B x T): 620 x 290 x 420mm - Weight: 68kg						