

xrPress ET25 - MOTORIZED Hydraulic Presses

Art. ID LM24-PRESS-ET25
Unit 68kg

Description

The MotorT25, is a power assisted hydraulic press operating to 25 Tons. The press has been designed to handle a wide variety of pressing applications, including XRF and IR sample preparation.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|---|-------------------|------------|---------------------|------|--------|--------|
| Key features | | | | | | |
| Large working distance between pressing surfaces | | | | | | |
| Low noise operation | | | | | | |
| Integral high clarity PETG safety Guards | | | | | | |
| Max. Load on Ram (Piston) n) T25: 25 Tons | | | | | | |
| Load Range (0.5Ton increments) T25: 2 – 25 Ton | | | | | | |
| Top Bolster Diameter: 400mm | | | | | | |
| Top Lead Screw Vertical travel: 90mm | | | | | | |
| Ram (Piston) Bolster Diameter: 100mm | | | | | | |
| Ram (Piston) Stroke: 25mm | | | | | | |
| Max/Min Dist. between Pressing Faces: 165mm – 70mm | | | | | | |
| Technical specifications | | | | | | |
| Dimensions (H x B x T): 620 x 290 x 420mm - Weight: 68kg | | | | | | |