

Custom Solution - 6 Components in 5% HNO₃,

Art. ID LM24-CUS-70508
Unit 3 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	1	mg/L		
	Boron (B)	[7440-42-8]	1	mg/L		
	Potassium (K)	[7440-09-7]	2	mg/L		
	Magnesium (Mg)	[7439-95-4]	10	mg/L		
	Sodium (Na)	[7440-23-5]	10	mg/L		
	Calcium (Ca)	[7440-70-2]	50	mg/L		