

Custom Solution - 17 Components in 5% HNO₃

Art. ID	LM24-CUS-70509
Unit	2 x 100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cadmium (Cd)	[7440-43-9]	5	µg/L		
	Silver (Ag)	[7440-22-4]	10	µg/L		
	Antimony (Sb)	[7440-36-0]	5	µg/L		
	Arsenic (As)	[7440-38-2]	10	µg/L		
	Nickel (Ni)	[7440-02-0]	20	µg/L		
	Lead (Pb)	[7439-92-1]	10	µg/L		
	Selenium (Se)	[7782-49-2]	10	µg/L		
	Chromium (Cr)	[7440-47-3]	50	µg/L		
	Barium (Ba)	[7440-39-3]	20	µg/L		
	Vanadium (V)	[7440-62-2]	50	µg/L		
	Manganese (Mn)	[7439-96-5]	50	µg/L		
	Aluminium (Al)	[7429-90-5]	200	µg/L		
	Iron (Fe)	[7439-89-6]	200	µg/L		
	Strontium (Sr)	[7440-24-6]	400	µg/L		
	Uranium (U)	[7440-61-1]	20	µg/L		
	Beryllium (Be)	[7440-41-7]	4	µg/L		
	Cobalt (Co)	[7440-48-4]	50	µg/L		