

PAK - Mix 7 - 100 µg/mL of each analyte in Methanol

Art. ID LM24-CUS-10123
Unit 6 x 2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fluoranthene	[206-44-0]	100	µg/mL		
	Benzo(b)fluoranthene	[205-99-2]	100	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	100	µg/mL		
	Benzo(a)pyrene	[50-32-8]	100	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	100	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	100	µg/mL		
	Tetrabutyltin	[1461-25-2]	100	µg/mL		