

## **High-quality nickel alloy, Chips**

Art. ID VS-N2/3  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,2	%		
	Arsenic (As)	[7440-38-2]	<0,001	%		
	Carbon (C)	[7440-44-0]	0,018	%		
	Chromium (Cr)	[7440-47-3]	5,59	%		
	Copper (Cu)	[7440-50-8]	0,083	%		
	Iron (Fe)	[7439-89-6]	Rem	%		
	Manganese (Mn)	[7439-96-5]	0,84	%		
	Nickel (Ni)	[7440-02-0]	76,3	%		
	Phosphorus (P)	[7723-14-0]	0,0034	%		
	Sulfur (S)	[7704-34-9]	0,0025	%		
	Silicon (Si)	[7440-21-3]	1,4	%		