

## High-quality nickel alloy, Chips

Art. ID VS-N4/3  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	<0,001	%		
	Carbon (C)	[7440-44-0]	0,0057	%		
	Chromium (Cr)	[7440-47-3]	0,07	%		
	Copper (Cu)	[7440-50-8]	5,65	%		
	Iron (Fe)	[7439-89-6]	5,8	%		
	Manganese (Mn)	[7439-96-5]	0,762	%		
	Molybdenum (Mo)	[7439-98-7]	4,87	%		
	Phosphorus (P)	[7723-14-0]	0,0019	%		
	Silicon (Si)	[7440-21-3]	0,81	%		