

**Cast iron and pig iron, truncated cone h 38-42mm, Ø of top base – 36-40mm, of bottom base – 38-42mm, Disc**

Art. ID VS-ChG10/6  
Unit Disc  
Deliverydetails No Dangerous Good

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	(0,003-0,006)	%		
	Carbon (C)	[7440-44-0]	~2,7	%		
	Chromium (Cr)	[7440-47-3]	~0,2	%		
	Manganese (Mn)	[7439-96-5]	0,86	%		
	Phosphorus (P)	[7723-14-0]	0,103	%		
	Sulfur (S)	[7704-34-9]	0,0072	%		
	Silicon (Si)	[7440-21-3]	2,86	%		
	Vanadium (V)	[7440-62-2]	~0,3	%		