

**Alloy cast iron, truncated tetrahedral pyramid h 20-26mm, side of working surface – 35-40 mm, Block**

Art. ID VS-ChG24  
Unit Block  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,007	%		
	Carbon (C)	[7440-44-0]	3,05	%		
	Chromium (Cr)	[7440-47-3]	0,031	%		
	Copper (Cu)	[7440-50-8]	0,1	%		
	Magnesium (Mg)	[7439-95-4]	0,015	%		
	Manganese (Mn)	[7439-96-5]	0,245	%		
	Molybdenum (Mo)	[7439-98-7]	0,031	%		
	Nickel (Ni)	[7440-02-0]	0,87	%		
	Phosphorus (P)	[7723-14-0]	0,26	%		
	Sulfur (S)	[7704-34-9]	0,0048	%		
	Antimony (Sb)	[7440-36-0]	0,009	%		
	Silicon (Si)	[7440-21-3]	2,5	%		
	Tin (Sn)	[7440-31-5]	0,077	%		
	Titanium (Ti)	[7440-32-6]	0,06	%		
	Vanadium (V)	[7440-62-2]	0,0067	%		