

Alloy cast iron, truncated cone h 35-40mm, Ø of top base – 32-36mm, of bottom base – 35-40mm, Disc

Art. ID VS-ChG35
Unit Disc
Deliverydetails No Dangerous Good

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,34	%		
	Chromium (Cr)	[7440-47-3]	0,233	%		
	Copper (Cu)	[7440-50-8]	0,09	%		
	Manganese (Mn)	[7439-96-5]	1,23	%		
	Molybdenum (Mo)	[7439-98-7]	0,027	%		
	Nickel (Ni)	[7440-02-0]	2,15	%		
	Phosphorus (P)	[7723-14-0]	0,102	%		
	Sulfur (S)	[7704-34-9]	0,021	%		
	Silicon (Si)	[7440-21-3]	0,617	%		
	Titanium (Ti)	[7440-32-6]	0,022	%		
	Vanadium (V)	[7440-62-2]	0,043	%		