

**Alloy cast iron, truncated cone h 35-40mm, Ø of top base – 32-36mm, of bottom base – 35-40mm, Disc**

Art. ID VS-ChG38  
Unit Disc  
Deliverydetails No Dangerous Good

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,43	%		
	Chromium (Cr)	[7440-47-3]	1,98	%		
	Copper (Cu)	[7440-50-8]	1,2	%		
	Manganese (Mn)	[7439-96-5]	0,302	%		
	Molybdenum (Mo)	[7439-98-7]	0,046	%		
	Nickel (Ni)	[7440-02-0]	0,162	%		
	Phosphorus (P)	[7723-14-0]	0,386	%		
	Sulfur (S)	[7704-34-9]	0,084	%		
	Silicon (Si)	[7440-21-3]	2,3	%		
	Titanium (Ti)	[7440-32-6]	0,105	%		
	Vanadium (V)	[7440-62-2]	0,119	%		