

Alloy cast iron, truncated cone h 35-40mm, Ø of top base – 32-36mm, of bottom base – 35-40mm, Disc

Art. ID VS-ChG40
Unit Disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,59	%		
	Chromium (Cr)	[7440-47-3]	1,47	%		
	Copper (Cu)	[7440-50-8]	0,98	%		
	Manganese (Mn)	[7439-96-5]	1,56	%		
	Molybdenum (Mo)	[7439-98-7]	0,229	%		
	Nickel (Ni)	[7440-02-0]	1,61	%		
	Phosphorus (P)	[7723-14-0]	0,059	%		
	Sulfur (S)	[7704-34-9]	0,019	%		
	Silicon (Si)	[7440-21-3]	1,6	%		
	Titanium (Ti)	[7440-32-6]	0,18	%		
	Vanadium (V)	[7440-62-2]	0,325	%		