

**Alloy cast iron, truncated tetrahedral pyramid h 20-25mm, Ø of top base – 28-34mm, of bottom base – 35-39mm, Block**

Art. ID VS-ChG47  
Unit Block  
Deliverydetails No Dangerous Good

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0056	%		
	Arsenic (As)	[7440-38-2]	0,014	%		
	Carbon (C)	[7440-44-0]	2,43	%		
	Cobalt (Co)	[7440-48-4]	0,0042	%		
	Chromium (Cr)	[7440-47-3]	14,45	%		
	Copper (Cu)	[7440-50-8]	0,0104	%		
	Manganese (Mn)	[7439-96-5]	0,949	%		
	Molybdenum (Mo)	[7439-98-7]	0,0019	%		
	Nickel (Ni)	[7440-02-0]	0,149	%		
	Phosphorus (P)	[7723-14-0]	0,099	%		
	Sulfur (S)	[7704-34-9]	0,083	%		
	Antimony (Sb)	[7440-36-0]	0,04	%		
	Silicon (Si)	[7440-21-3]	2,73	%		
	Tin (Sn)	[7440-31-5]	0,093	%		
	Titanium (Ti)	[7440-32-6]	0,041	%		
	Vanadium (V)	[7440-62-2]	0,129	%		