

Custom Solution - Chromium (Cr) - 100 mg/L in HNO₃ for ICP-MS

Art. ID LM24-CUS-70522
Unit 1 x 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle) / (Cr(NO₃)₃ 99.999%)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	100	mg/L		