

### **Custom Solution - Lead (Pb) - 100 mg/L in HNO3 for ICP-MS**

Art. ID LM24-CUS-70523  
Unit 1 x 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation: 24 months (originally closed bottle) / (Pb(NO<sub>3</sub>)<sub>2</sub> 99.9995%).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	100	mg/L		