

Custom Solution - Nickel (Ni) - 100 mg/l in HNO3 for ICP-MS

Art. ID LM24-CUS-70524
Unit 1 x 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle) / (Ni, 99.9998%)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nickel (Ni)	[7440-02-0]	100	mg/L		