

Custom Solution - Arsenic (As) - 100 mg/L in HNO3 for ICP-MS (As 99.999%).

Art. ID LM24-CUS-70525
Unit 1 x 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle) / (As 99.999%)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		