

### **Custom Solution - 14 components in HNO3 5%**

Art. ID LM24-CUS-70534  
Unit 1 x 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation: 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	10	mg/L		
	Arsenic (As)	[7440-38-2]	25	mg/L		
	Barium (Ba)	[7440-39-3]	500	mg/L		
	Boron (B)	[7440-42-8]	1000	mg/L		
	Cadmium (Cd)	[7440-43-9]	10	mg/L		
	Lithium (Li)	[7439-93-2]	50	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Lead (Pb)	[7439-92-1]	25	mg/L		
	Copper (Cu)	[7440-50-8]	250	mg/L		
	Selenium (Se)	[7782-49-2]	5	mg/L		
	Thallium (Tl)	[7440-28-0]	5	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	50	mg/L		