

### **Custom Solution - CKW Mix 2 je 2000 µg/mL in Methanol**

Art. ID LM24-CUS-10126  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,1,2-Trichloroethane	[79-00-5]	2000	µg/mL		
	1,1,2,2-Tetrachloroethane	[79-34-5]	2000	µg/mL		