

Carbon steel, 0.2 C (chips)

Art. ID	NIST-19h
Unit	150 g
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is in the form of chips sized between 0.50 and 1.18 mm sieve openings (35 and 16 mesh). Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,215	%		
	Manganese (Mn)	[7439-96-5]	0,393	%		
	Phosphorus (P)	[7723-14-0]	0,016	%		
	Sulfur (S)	[7704-34-9]	0,022	%		
	Silicon (Si)	[7440-21-3]	0,211	%		
	Copper (Cu)	[7440-50-8]	0,466	%		
	Nickel (Ni)	[7440-02-0]	0,248	%		
	Chromium (Cr)	[7440-47-3]	0,173	%		
	Vanadium (V)	[7440-62-2]	0,003	%		
	Molybdenum (Mo)	[7439-98-7]	0,038	%		
	Aluminium (Al)	[7429-90-5]	0,002	%		