

2-Amino-5-nitroso-4,6-pyrimidinediol

Art. ID TRC-A618735-500MG
Unit 500 mg
Deliverydetails /not restricted

Description

Application Note: 2-Amino-5-nitroso-4,6-pyrimidinediol is used as a reagent in the synthesis of 6-amino-2,3-dimethylpyrimido[4,5-e][1,2,4]triazine-8-one as an inhibitor of dihydropteridine reductase (DHPR). 2-Amino-5-nitroso-4,6-pyrimidinediol is also used in analytical studies for determination of ruthenium and transition metals.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Amino-5-nitroso-4,6-pyrimidinediol	[55482-22-9]				