

(S)-3-Amino(piperidine-3-d1)-2,6-dione hydrochloride

Art. ID TRC-A609132-10MG
Unit 10 mg

Description

Application Note: (S)-3-Amino(piperidine-3-d1)-2,6-dione HCl is deuterium labelled (S)-3-Amino(piperidine-3)-2,6-dione HCl, which can be used to prepare thalidomide analogs and lenalidomide.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(S)-3-Amino(piperidine-3-d1)-2,6-dione HCl	[1398112-31-6]				