

N-Acetyl-S-phenyl-DL-cysteine-3,3-d2

Art. ID TRC-A187359-2.5MG

Unit 2,5 mg

Description

Application Note: N-Acetyl-S-phenyl-DL-cysteine-3,3-d2 is the deuterium labelled form of N-Acetyl-S-phenyl-DL-cysteine, which is used as a reactant in the preparation of potent HIV-1 protease inhibitors that use D-amino acids as P2/P3 ligands.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-S-phenyl-DL-cysteine-3,3-d2					