

## **Ferrosilicon**

Art. ID                      NIST-59a  
Unit                         50 g  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Cargo Air Only (CAO) UN1408 (6.1) III E1

### Description

Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	48,1	%		
	Iron (Fe)	[7439-89-6]	50,05	%		
	Carbon (C)	[7440-44-0]	0,046	%		
	Manganese (Mn)	[7439-96-5]	0,75	%		
	Phosphorus (P)	[7723-14-0]	0,016	%		
	Sulfur (S)	[7704-34-9]	0,002	%		
	Copper (Cu)	[7440-50-8]	0,052	%		
	Nickel (Ni)	[7440-02-0]	0,033	%		
	Chromium (Cr)	[7440-47-3]	0,08	%		
	Aluminium (Al)	[7429-90-5]	0,35	%		
	Boron (B)	[7440-42-8]	0,058	%		
	Calcium (Ca)	[7440-70-2]	0,042	%		