

**Iron base, Steel - Stainless, CRM, Disc 38 mm x 30 mm, please see NCS HS11709A-1--7**

Art. ID                      NCS HS11709a-6  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 0,174               | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,0066              | %    |        |        |
|                  | Carbon (C)        | [7440-44-0] | 0,027               | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 0,031               | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 24,4                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,326               | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,595               | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 0,053               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 7,34                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,057               | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,0002              | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,042               | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 1,06                | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,042               | %    |        |        |
|                  | Titanium (Ti)     | [7440-32-6] | 0,17                | %    |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 0,078               | %    |        |        |